


FILLING YOUR NEEDS



BOTTELPOST

UNIT DOSES


FOR OPHTHALMIC APPLICATIONS/PRODUCTS



bottelpack®

BLOW-FILL-SEAL TECHNOLOGY


Detach one unit
from the card.



Hold it and turn
the head flag off.




Apply the product
by squeezing the ampoule.



Discard any unused product.





Ampoules can be presented as unit doses and multi doses. They are safe for administration and are used for a broad spectrum of products, such as injectables, inhalation therapy, droppers, irrigation and other.

For the double administration of eye drops, nose drops and ear drops, unit dose ampoules also with reclosable twist-off design have been developed. There is no need for preservatives.

The advantages of unit doses are obvious:

- No preservatives - no allergic reactions
- No cross contamination - one user only
- Easy handling with correct dosing
- Environmentally friendly packing material

The Twist-off Closure is utilized by the millions as a very economical closure which can be used with many products.

ROMMELAG ENGINEERING

Rommelag Kunststoff-Maschinen
Vertriebsgesellschaft mbH

Mayenner Straße 18-20
71332 Waiblingen
Germany

T + 49 7151 95811-0

mail.rd@rommelag.com
www.rommelag.com

Rommelag Trading (Shanghai) Co., Ltd.

Room 104~107, 1st Floor · No. 510
Zhongshan South 2nd Road · Xuhui District
200032 Shanghai · P.R. China

T + 86 21 515 349 67

mail.rcn@rommelag.com
www.rommelag.com

Rommelag AG

Fabrikweg 16
5033 Buchs
Switzerland

T + 41 62 83455-55

mail.rch@rommelag.com
www.rommelag.com

Rommelag USA, Inc.

27905 Meadow Drive
Suite 9 · Evergreen
CO 80439 · USA

T + 1 303 674 8333

mail.romus@rommelag.com
www.rommelag.com